

# TRIAL LENS SET SERIES



### Specifications

TRIAL LENS SET TL-34P (PLASTIC RIM) 232PCS									
CONCAVE(±) SPHERE				CONCAVE(±) CYLINDER		PRISM		ACCESSORY	
LIST	PCS	LIST	PCS	LIST	PCS	LIST	PCS	SPECIFICATION	PCS
0.12	2	5.00	2	0.12	2	1.00	1	Maddox	1
0.25	2	5.50	2	0.25	2	2.00	1	Frosted	1
0.50	2	6.00	2	0.50	2	3.00	1	Pin Hole (1mm)	1
0.75	2	6.50	2	0.75	2	4.00	1	Occulder	1
1.00	2	7.00	2	1.00	2	5.00	1	Slit	1
1.25	2	8.00	2	1.25	2	6.00	1	Red Filter	1
1.50	2	9.00	2	1.50	2	7.00	1	Green Filter	1
1.75	2	10.00	2	1.75	2	8.00	1	Plane	1
2.00	2	11.00	2	2.00	2	9.00	1	Crossed Line	2
2.25	2	12.00	2	2.25	2	10.00	1		
2.50	2	13.00	2	2.50	2				
2.75	2	14.00	2	2.75	2				
3.00	2	15.00	2	3.00	2				
3.25	2	16.00	2	3.25	2				
3.50	2	18.00	2	3.50	2				
3.75	2	20.00	2	4.00	2				
4.00	2			4.50	2				
4.50	2			5.00	2				
				6.00	2				

TRIAL LENS SET TL-35M (METAL RIM) 226PCS									
CONCAVE(±) SPHERE				CONCAVE(±) CYLINDER		PRISM		ACCESSORY	
LIST	PCS	LIST	PCS	LIST	PCS	LIST	PCS	SPECIFICATION	PCS
0.12	2	5.00	2	0.12	2	0.50	2	Pin Hole (1mm)	2
0.25	2	5.50	2	0.25	2	1.00	2	Occulder	1
0.50	2	6.00	2	0.50	2	2.00	2	Slit	1
0.75	2	6.50	2	0.75	2	3.00	1	Red Filter	1
1.00	2	7.00	2	1.00	2	4.00	1	Green Filter	1
1.25	2	7.50	2	1.25	2	5.00	1		
1.50	2	8.00	2	1.50	2	6.00	1		
1.75	2	9.00	2	1.75	2	7.00	1		
2.00	2	10.00	2	2.00	2	8.00	1		
2.25	2	11.00	2	2.25	2				
2.50	2	12.00	2	2.50	2				
2.75	2	13.00	2	3.00	2				
3.00	2	14.00	2	3.50	2				
3.25	2	15.00	2	4.00	2				
3.50	2	16.00	2	4.50	2				
3.75	2	18.00	2	5.00	2				
4.00	2	20.0	2	6.00	2				
4.50	2								

TRIAL LENS SET TL-24P(PLASTIC RIM) 152PCS									
CONCAVE(±) SPHERE				CONCAVE(±) CYLINDER		PRISM		ACCESSORY	
LIST	PCS	LIST	PCS	LIST	PCS	LIST	PCS	SPECIFICATION	PCS
0.25	2	3.50	2	0.25	2	0.50	2	Occulder	1
0.50	2	3.75	2	0.50	2	1.00	2	Pin Hole (1mm)	1
0.75	2	4.00	2	0.75	2	2.00	1	Slit	1
1.00	2	4.50	2	1.00	2	3.00	1	Maddox	1
1.25	2	5.00	2	1.25	2	4.00	1	Crossed Line	2
1.50	2	5.50	2	1.50	2	5.00	1	Red Filter	1
1.75	2	6.00	2	1.75	2			Green Filter	1
2.00	2	7.00	2	2.00	2				
2.25	2	8.00	2	2.50	2				
2.50	2	9.00	2	3.00	2				
2.75	2	10.00	2						
3.00	2								
3.25	2								

TRIAL LENS SET TL-24M (METAL RIM) 132PCS									
CONCAVE(±) SPHERE				CONCAVE(±) CYLINDER		PRISM		ACCESSORY	
LIST	PCS	LIST	PCS	LIST	PCS	LIST	PCS	SPECIFICATION	PCS
0.12	2	3.00	2	0.25	2			Occulder	1
0.25	2	3.25	2	0.50	2			Pin Hole (1mm)	1
0.50	2	3.50	2	0.75	2			Slit	1
0.75	2	3.75	2	1.00	2			Plane	1
1.00	2	4.00	2	1.25	2				
1.25	2	4.50	2	1.50	2				
1.50	2	5.00	2	1.75	2				
1.75	2	5.50	2	2.00	2				
2.00	2	6.00	2						
2.25	2	7.00	2						
2.50	2	8.00	2						
2.75	2	10.00	2						

TRIAL LENS SET TL-28M 68PCS					
CONCAVE(±) SPHERE				CONCAVE(-) CYLINDER	
LIST	PCS	LIST	PCS	LIST	PCS
0.25	1	4.50	1	0.25	1
0.50	1	5.00	1	0.50	1
0.75	1	5.50	1	0.75	1
1.00	1	6.00	1	1.00	1
1.25	1	6.50	1	1.50	1
1.50	1	7.00	1	2.00	1
1.75	1	8.00	1	2.50	1
2.00	1	9.00	1	3.00	1
2.25	1	10.00	1	3.50	1
2.50	1	11.00	1	4.00	1
2.75	1	12.00	1	4.50	1
3.00	1	13.00	1	5.00	1
3.50	1	14.00	1		
4.00	1	16.00	1		

### Trial Frame



**TF-1/TF-1A**  
 PD: 48mm - 80mm  
 NET WEIGHT : 85g (TF-1)  
 NET WEIGHT : 90g (TF-1A)



**TF-3**  
 PD: 48mm - 80mm  
 NET WEIGHT : 60g



**Light Weight Trial Frame**  
 PD: 56, 58, 60, 62, 64, 66mm  
 NET WEIGHT : 39.5g - 40.5g

Design and specifications are subject to change without notice.

Manufacturer

**Rexxam**  
 Rexxam Co.,Ltd.

Kagawa factory  
 958, Ikeuchi, Konan-cho,  
 Takamatsu, Kagawa 761-1494, Japan

Contact

**Eye-care Instruments Sales Dept. Tokyo Office**

2-4-2 Kandatsukasa-machi, Chiyoda-ku, Tokyo, 101-0048, Japan

TEL:81-3-6262-9471 FAX:81-3-6262-9472

E-mail: eye@rexam.co.jp URL : http://www.rexxam.co.jp URL : http://www.shin-nippon.co.jp

Distributed by



Management System  
 ISO 9001:2008  
 ISO 13485:2003

www.tuv.com  
 ID 0910589071

